

# Deluxe rechargeable HANDHELD PULSE OXIMETER

When you are looking for a deluxe handheld oximeter at an affordable price, you don't want to give up or compromise on quality. This deluxe rechargeable handheld pulse oximeter gives accurate readings in a durable unit that you can rely on, patient after patient. Convenient waveform display option helps you better assess the patient's condition. This convenient, portable monitor is ideal for spot checks or continuous monitoring in virtually any setting.

- Designed with curvature that fits easily into the nurses hand.
- Lightweight and easy to use
- Built-in rechargeable battery and power save function
- 1.8" high definition LCD color screen with a realtime showing of SpO2 waveform, SpO2 & PR value
- Unique digital SpO2 technology provides fast, accurate, stable, and strong anti-interference performance
- Real time clock with easy setup
- **BCI compatible finger probe**
- 10-12 hours battery back-up with a real time battery charge indicator
- Built in memory card for data memory
- Suitable for adults, pediatric, and neonates.
- Suitable for use in nursing homes, physicians offices, clinics, outpatient surgical centers, long term care facilities, & any other setting looking for reliable & affordable patient monitoring.
- Able to do spot or continuous monitoring



## Technical Specifications

### SpO2

**Patient type:** Adult, pediatric, neonate  
**Measurement range:** 0-100%  
**Resolution:** 1%  
**Accuracy:**  $\pm$  2bpm  
**Parameter:** SpO2, PR

### Pulse Rate(PR)

**Measurement range:** 30-250bpm  
**Resolution:** 1bpm  
**Accuracy:**  $\pm$  2bpm

### Work environment

**Temperature range operating:** 0°C - 50°C  
**Transport/Storage:** -40°C - 131°C

### Relative Humidity

**Operating:** 0%-93%  
**Transport/Storage:** 0%-95%

### Power requirements

**Operating:** 100-240VAC, 50/60HZ  
**Power consumption:** <6VA



**ITEM # 20120**

*Digital and  
Waveform  
Display Options*



 **Proactive**  
medical products

**270 Washington Street, Mount Vernon, NY 10553**  
**Tel: 855-BE-PROACTIVE (855-237-7622) Fax: 845-205-6005**  
**www.proactivemedical.com • sales@proactivemedical.com**