

## DermaFloat® LAL Therapeutic Surface

Technology that will revolutionize your wound prevention and intervention

- True low air loss wicks away moisture to prevent skin maceration
- Twenty individual cells help to evenly distribute resident's weight and maximize pressure redistribution
- Low shear and low friction therapy pad is treated with an antibacterial and antifungal agent and is fluid proof and moisture vapor permeable
- Two-inch lower safety mattress for resident transport and power failure; stays inflated for up to 12 hours to prevent bottoming out
- CairRails risk management side air bolsters



The DermaFloat® LAL, provided by Joerns Healthcare is a unique therapy system that provides pressure redistribution by combining low air loss with pulsation.

Low air loss therapy has been demonstrated to reduce the risk of pressure ulcers as well as being a valuable aid in the treatment of pressure ulcers. The pulsation feature may provide additional resident comfort and enhance the lymphatic system. The DermaFloat LAL is designed as a mattress replacement system for residents weighing up to 500 lbs.



Eight comfort control settings to maximize resident compliance and promote healing

## Technical Specifications

Feature	
Maximum Weight Capacity	500 lbs.
Dimensions	36" W x 80" L x 10" D 42" W x 80" L x 10" D
Number of Zones	3
Number of Therapy Cells	20 horizontal cells
Depth of Therapy Cells	8"
Lower Safety Cell Depth	2"
CPR Quick Deflation	Yes
CairRails/Side Bolsters	CairRails
Control Unit Dimensions	7.5"W x 12.25"H x 5.5"D
Control Unit Modes of Operation	Autofirm, Therapy and Pulsate
Therapy / Comfort Settings	8
Standby	Yes
Audible Power Failure Alarm	Yes
Visual Power Failure Alarm	Yes
UL / CUL Classified / Listed	UL 60601-1 CAN/CSA C22.2 No. 60601-1
Warranty	2 years: mattress 1 year: cover 1 year: electromechanical mattress components

Joerns is revolutionizing healthcare through design, prevention and advanced healing solutions. We are focused on improving patient and caregiver experiences, placing the well being of others at the center of everything we do. To maximize the performance of our customers, we bring to market fully integrated services and technologies that provide the highest quality and represent the greatest value. For more information visit [www.joerns.com](http://www.joerns.com).