

# Joerns BariMatt 1000 Plus Therapeutic Support Surface



## BariMatt 1000 Plus

The BariMatt 1000 Plus is a uniquely designed expandable foam bariatric mattress that provides caregivers six different configurations to meet the needs of all of their bariatric patients. This flexibility can help reduce costs while ensuring optimal clinical outcomes.

Designed to enhance the comfort of the Joerns Bari10A bed frame, the convenient side expansion panels conveniently zip on and off to allow ease of patient care and transport.

BariMatt 1000 Plus mattress and bed compatibility 36"W x 80"L x 6.5"D 42"W x 80"L x 6.5"D



Shown in smallest configuration.



Shown with 6" expansion panel inserted.



48"W x 88"L x 6.5"D

Shown with two 6" expansion panels inserted and 8" length extender.

#### **Addressing Bariatric Challenges**

With 133.6 million (66.3%) of Americans overweight and 32% of those being clinically obese, many facilities are being forced to redesign their existing room environments to accommodate these patients. This expense poses great financial challenges for our healthcare systems.<sup>4</sup>

Joerns Healthcare can help facilities to address these rising costs and provide facilities with solutions to address their bariatric needs by providing therapeutic support surfaces and bariatric bed frames that can:

- Improve patient comfort, healing and mobility
- Help treat and prevent pressure ulcers
- · Manage risks associated with bariatric care to both the patient and caregiver

<sup>1</sup>Center for Disease Control (http://www.cdc.gov/nchs/fastats/overwt.htm)

# transforming wound management


### **Key features**

All foam elements are enclosed with durable, moisture resistant covers for easy cleaning and to protect the foam when not installed.

- Available in three width configurations: 36", 42" and 48" and two length configurations: 80" and 88"
- Convenient transport mode (36" W x 80" L) allows the mattress size to be adjusted to easily navigate through standard doorways and elevators while the patient is on the mattress
- · Mattress relief cuts allowing the surface to articulate with the bed frame for enhanced patient comfort
- Fabric mattress surround is designed to neatly fit together and without gaps providing patients with a comfortable, safe healing environment
- · One-inch visco-elastic foam topper envelops the patient for optimized comfort
- · Top cover has screened illustrations allowing caregivers to quickly set up the mattress and meet patient's individual needs

Removable foot and side bolsters





Tied support straps





Screened mattress instructions



## Joerns Bari10A Bed

The Bari10A, as it is commonly called, represents a revolutionary new concept in bariatric beds. The newly designed Bari10A was specifically developed to address patient dignity in hospital and rehabilitation facilities while providing safety and comfort when weight exceeds standard surface weight limits.

In addition, the Bari10A has a fully-articulated mattress support platform with integrated width and length extension for comfort and transport modes. The Bari10A bed design allows for independent head and knee movement, up/down adjustments, and Trendelenberg and Reverse Trendelenberg positioning, while ensuring optimal clinical outcomes for bariatric patients.

"With decades of experience, Joerns Healthcare understands the healthcare bed market. We have developed the new Joerns Bari10A to complete your healthcare bed needs," said J.B. Risk, Joerns Senior Product Manager. "It is designed for a 1050 lb. work load and is sturdy throughout, with four motors providing independent functions. We are confident in the design, the production and strength of our new bariatric bed product."

Joerns Healthcare is a leading manufacturer and trusted source of advanced products and services to the healthcare market. Industry-leading brands include Joerns beds and furnishings, TriLine and BioClinic wound management and Hoyer® lifting and repositioning products and services. Additionally, value-added programs are available as a part of EquipTech Solutions, Creating a Safe Environment, and Woundlink.





Post Acute, Acute, HomeCare 2430 Whitehall Park Drive Charlotte, NC 282731 (P) 800.826.0270 (F) 800.457.8827 VA/Government 19748 Dearborn St Chatsworth, CA 91311 (P) 800.966.6662 (F) 800.232.9796 Canadian Office Clarke Road, Suite 6 London, ON Canada N5V 3A9 (P) 866.546.1151 (F) 519.451.8662

© 2012 Joerns Healthcare • 6220046 RevB • 12-0055