



Protekt™ Gel Overlay



 **Proactive**
medical products

Medicare
Code
EO185

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com

The Protekt™ Gel Overlay utilizes a series of gel bladders encapsulated between high density, high resilient, contoured foam and a firm support base. This composition creates a cost effective pressure redistribution group I support surface.

FEATURES AND BENEFITS:

- Optimal overlay surface for the prevention and treatment of pressure ulcers.
- Convulated high density polyurethane foam top layer.
- Bladders contain an aqueous gel composition that provides deep submersion.
- Base foundation is a high density high resilient foam construction which prevents “bottoming out”.
- Cover is a soft nylon urethane coated vapor permeable water proof fabric.
- Top cover is non-shear, low-friction, and anti-microbial/bacterial.
- 4 elastic corner straps secures overlay to mattress.
- Gel maintains room temperature at all times.
- Bariatric sizes available.
- 16 CFR part 1632 & 1633.
- Weight capacity: 300 lbs.
- Size: 34"x78"x3.5".
- Warranty: one year non-prorated.

90001	Gel Overlay	Gel Overlay (EO185)	34"x78"x3"
90002	Gel Overlay 4200	Gel Overlay	42"x78"x4"
90003	Gel Overlay 4200B	Gel Overlay Bariatric	48"x78"x4"
90004	Gel Overlay 4800B	Gel Overlay Bariatric	54"x78"x4"



270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com