



SC-3004-DL Sequential Circulator

Model SC-3004-DL Sequential Circulator is a gradient, sequential, pneumatic compression device, intended for the primary or adjunctive treatment of primary or secondary Lymphedema. The device is additionally intended for alternate treatment of chronic venous stasis ulcers and associated venous insufficiency, as well as general treatment of swelling of the extremities. The SC-3004-DL is intended for home or hospital use.

Key Features:

- Unilateral and/or Bilateral Functionality
- Non-peristaltic
- Custom Garment Availability
- Light Weight and Portable
- Easy to Use
- 3 Year Warranty

Benefits:

- Reduces Pain
- Improves Lymphatic and Vascular Flow
- Promotes Wound Healing
- Prevents Infection
- Provides Effective Prophylaxis for Venous Thrombosis
- Reduces Dynamic Edema



SUPPLY ITEM: E0652



Contraindications:

Pneumatic compression is contraindicated for patients with the following:

- Congestive heart failure
- Deep Vein Thrombosis
- Inflammatory phlebitis or episodes of pulmonary embolism
- Infections in the limb, including cellulitis, without appropriate antibiotic coverage
- Active Cancer except for palliative care
- Other indications as identified by the treating physician

WARNING: Federal law restricts this device to sale by or on the order of a physician.



Bio Compression Systems' Biocomfort Garments

Made of high quality 200 Denier Nylon Oxford coated with 3 mils of Polyurethane
Available in a variety of sizes and styles including: over the shoulder,
adjustable, custom and amputee.

| Upper Extremity/4-Chamber Garment Sizing Chart | | | | | | | | |
|--|--------|----------------|-------|-------|---------|---------|----------------------|------------------------|
| Part Number | Size | Circumferences | | | | Lengths | | |
| | | Axilla | Bicep | Elbow | Forearm | Wrist | Axilla to Fingertips | Shoulder to Fingertips |
| GS-3035-S | Small | 23.5" | 22" | 19.5" | 17.5" | 15" | 24" | 26.5" |
| GS-3035-M | Medium | 27" | 26" | 24" | 21" | 18" | 28" | 31" |
| GS-3035-L | Large | 25.5" | 24" | 21.5" | 19" | 16" | 31" | 35.5" |
| GS-3035-SH/S | Small | 23.5" | 22" | 19.5" | 17.5" | 15" | 24" | 26.5"+8" |
| GS-3035-SH/M | Medium | 27" | 26" | 24" | 21" | 18" | 28" | 31"+8" |
| GS-3035-SH/L | Large | 25.5" | 24" | 21.5" | 19" | 16" | 31" | 35.5"+8" |

NOTE: Measurements above are in **Inches** and are the maximum allotted point on limb

| Lower Extremity/4-Chamber Garment Sizing Chart | | | | | | | | | | | | |
|--|-----------------|--------------------|-----------|----------|----------|-------|---------------|-------------------------|-----------------|-----|-----|-----|
| Part Number | Size | Circumferences | | | | | Lengths | | | | | |
| | | Thigh below Crotch | Mid Thigh | Mid Knee | Mid Calf | Ankle | Arch (Instep) | Crotch to Heel (Length) | Inguinal Crease | | | |
| GS-3045-H | Half | Standard | N/A | 29" | 27" | 23.5" | 22" | 19" | NA | | | |
| GW-3045-H | | Wide | | 30.5" | 28" | 29" | 23" | | | | | |
| GA-3045-H | | Adjustable | | 32" | 30" | 23.5" | 22" | | | | | |
| GWA-3045-H | | Adj. Wide | | 33.5" | 31" | 29" | 23" | | | | | |
| GN-3045-S | Small | Narrow | 25.5" | 25" | 22" | 20" | 18" | 19.5" | 25.5" | 27" | | |
| GS-3045-S | | Standard | 31" | 30" | 28" | 24.5" | 21.5" | 19.5" | | | | |
| GW-3045-S | | Wide | 35.5" | 34" | 30.5" | 28" | 25.5" | 21.5" | | | | |
| GA-3045-S | | Adjustable | 34" | 33" | 31" | 27.5" | 21.5" | 19.5" | | | | |
| GWA-3045-S | | Adj. Wide | 38.5" | 37" | 33.5" | 31" | 25.5" | 21.5" | | | | |
| GN-3045-M | Medium | Narrow | 27.5" | 25.5" | 23" | 20" | 18" | 19.5" | 28" | 32" | | |
| GS-3045-M | | Standard | 32" | 30" | 28" | 24.5" | 21.5" | 19.5" | | | | |
| GW-3045-M | | Wide | 36.5" | 34" | 30" | 28" | 25.5" | 22.5" | | | | |
| GXW-3045 | | Extra Wide | 44" | 41" | 36" | 34" | 30" | 22.5" | | | 29" | 32" |
| GA-3045-M | | Adjustable | 35" | 34" | 31" | 27.5" | 21.5" | 19.5" | | | 28" | 32" |
| GWA-3045-M | | Adj. Wide | 39.5" | 36.5" | 33.5" | 31" | 25.5" | 22.5" | | | 28" | 32" |
| GXWA-3045-M | Adj. Extra Wide | 47" | 44" | 39" | 37" | 30" | 22.5" | 29" | 32" | | | |
| GN-3045-L | Large | Narrow | 28" | 26" | 24" | 21" | 18" | 19.5" | 32" | 36" | | |
| GS-3045-L | | Standard | 33" | 32" | 30" | 26" | 21.5" | 19.5" | | | | |
| GW-3045-L | | Wide | 37" | 36" | 32" | 29" | 24" | 20.5" | | | | |
| GA-3045-L | | Adjustable | 37" | 36" | 33" | 29" | 21.5" | 19.5" | | | | |
| GWA-3045-L | | Adj. Wide | 40" | 38" | 35" | 32" | 24" | 22.5" | | | | |

NOTE: Measurements above are in **Inches** and are the maximum allotted point on limb

