

PressureGuard® APM² Safety Supreme



Therapy Provided

- ✓ Powered Pressure Redistribution
- ✓ Alternating Pressure
- ✓ Lateral Rotation
- ✓ Raised Safety Perimeter

Combines all the features of the top-selling PressureGuard® APM², plus the added safety benefits of raised perimeter edges. The Safety Supreme incorporates the proprietary Geo-Mattress with Wings® raised perimeter design to help prompt the user away from the edge of the bed. Provides protection from pressure injuries, falls, fire, restraints and siderail entrapment – all in one easy-to-use package.

QUICK FACTS


Indications:

(in addition to APM² p. 8):

- Patients at risk for falls

Exclusive Raised Perimeter

incorporates the unique "edge-to-edge" design to help prompt user away from edge of bed without compromising skin integrity.



Item #	Description
SAF5875-29	Safety Supreme, 75"L x 35"W x 7"H center, 9"H side
SAF5880-29	Safety Supreme, 80"L x 35"W x 7"H center, 9"H side
SAF5884-29	Safety Supreme, 84"L x 35"W x 7"H center, 9"H side
5900	Digital control unit

Medicare Approved:	E0277 (Group II)
Weight limit:	350 lbs.
Warranty:	18 months (all components)
Mattress Weight:	approx. 20 lbs.
UL listed:	Yes

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.

