

PressureGuard® Bariatric APM

Therapy Provided

- ✓ Powered Pressure Redistribution
- ✓ Alternating Pressure for users up to 750 lbs.



QUICK FACTS

Indications:

- Aggressive treatment (through Stage 4)
- Prevention of pressure injuries in immobile or high risk patients
- Can't or won't reposition frequently
- Can be positioned off the pressure injury

Key Advantages:

- Simple operation
- Gentle and effective pressure changes
- Customizable firmness and time settings
- PressureGuard design: stability, safety, anti-shearing

Segmented top reduces sharp pressure drop-offs, while air system below handles the “heavy lifting,” alternating weight-bearing across the surface in user-selectable intervals. “Two up, one down” inflation pattern ensures maximum patient support throughout the cycle.

Item

Description

AP8042-29	Bariatric APM, 80”L x 42”W x 7”H
AP8048-29	Bariatric APM, 80”L x 48”W x 7”H
AP8054-29	Bariatric APM, 80”L x 54”W x 7”H
* AP8442-29	Bariatric APM, 84”L x 42”W x 7”H
AP8448-29	Bariatric APM, 84”L x 48”W x 7”H
* AP8454-29	Bariatric APM, 84”L x 54”W x 7”H
7500	Bariatric APM Control Unit Only

Medicare Approved: E0277 (Group II)

Weight limit: 750 lbs.

Warranty: 18 months (all components)

Mattress Weight: Approx. 25 lbs.

UL listed: Yes

Additional bariatric options available – see page 20.

All models are fire barrier version (FBV) as indicated by “-29” suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.

* Allow up to 10 working days lead time for this size.

Powerful Digital Control Unit

- Up to four times the air output of conventional units.
- Eight electronically-controlled comfort settings.
- Programmable cycle times of 20, 25, 30 or 35 minutes.
- Low pressure indicator and audible alarm.

