

OB-2600

ObboMed[®] Medical Air Mattress

Suitable for Stage 1-4 Pressure Ulcer

ObboMed[®] StandardAir[™] Mattress

10" Alternating with Low Air Loss • 15° Self-Lateral Rotation • 13" Side Bolster
Remote Control • 24-hour Power Outage Protection
(as long as air hose remains connected)

10" DEEP
2 LAYER MATTRESS

15° SELF-ROTATION

EASY TO OPERATE
CONTROL PANEL

REMOTE
CONTROL

13"H SIDE BOLSTER



ObboMed[®] is a global trade name and specializes in manufacturing the Medical air mattress; Air Pressure support surface mattress; and Alternating and Rotation air mattress. A professional team with 30 years of experience in medical devices—from design and production to marketing—We have been granted 15 invention patents covering USA, European Union, Japan, Russia, China, Taiwan, and pending in India.



obbomed.com

The ObboMed® difference- Comfort, Safety & Ease

Control Functions:

On/Off (Standby)
Lock
Max-Inflate (Auto Firm)
Static
Weight Display
Pressure Display
Alternating
Cycle Time
Power Alarm
Power Outage Protection
Low Air Loss
Low Air Pressure Warning
Fowler (Seating)

Additional Model Functions:

Side Bolster
Remote Control
15° Self-Lateral Rotation*



*Patented Features
For details, visit obbomed.com/patents

Features of ObboMed's multi-layered mattress:

- Designed to minimize friction and shearing of the skin
- Vapor permeable to facilitate moisture removal into the airflow away from the mattress
- Perspiration vacates through the mattress cover surface and/or the patient's bed boundary
- Antimicrobial, Liquid resistant, Fireproof
- Latex, lead, and mercury-free
- Below the cover layer, vapor flows out with the Low Air Loss Feature to lower humidity under the patient while simultaneously transferring away moisture and heat
- Air cells are composed of a combination of thermoplastic polyurethane and nylon (TPU) and airtight with high-frequency, heat sealing and are liquid proof and machine washable



**15° Self-Rotation By
Cone-Shaped Air Cells*
Low Air Loss up to 16 LPM**

OB-2600 Mattress dimension and weight specifications:
2-layer Air Cell, 10" Deep Mattress, 13"H Side Bolster

Mattress Dimension: 80"L x 36"W x 13"H
Pump: 20 LPM
Weight Capacity: 450 lbs
Control Unit Weight: 8.8 lbs

Mattress Weight: 24.5 lbs
Unit Package Weight: 37 lbs
Unit Package Dimension: 21"L x 18.5"W x 20.5"H

Limited Warranty:

1 year for Mattress and Top Cover, 2 years for Control Unit

©2022 ObboMed Group LLC. All rights reserved. Duplication, republication, or modification in whole or partially is prohibited without prior written permission from ObboMed. All trademarks are owned by or licensed to ObboMed Group LLC or its subsidiaries unless otherwise noted.



Basic Features

Therapy of Alternating with 15° Self-Rotation* and Low Air Loss for stages 1-4 pressure ulcer

24-hour Power Outage Protection (bottom layer stays inflated if air hose remains connected to Control Unit)

Fowler/Seating is designed to increase support in the seat/torso section of the mattress. When the head of the bed is raised, The fowler is intended to prevent any potential risk of patients from "bottoming out."

Sets of audio and visual warning signals if power outage and/or when air pressure is low

System defaults to "Alternating" mode automatically after Max inflate or 1-hour in Static Mode

User friendly control panel with Remote

Patented Features

Patented 15° Self-lateral Rotation is created automatically by cone-shaped air cell while in "alternating" mode.

PATENT INFORMATION:
<http://obbomed.com/patents>