



Innovative • Effective • Driven

Introducing our new line of Pneumatic Compression Devices.

With the same gradient sequential pneumatic compression we have always provided – now in a new, user-friendly digital version.

Use of the Bio Compression Systems pumps will not only manage your symptoms, it will improve your quality of life.



SC-2004-DL (E0651)

- 4 chamber compression therapy
- Unilateral or bilateral capability
- Adjustable pressure between 20 mmHg and 120 mmHg
- Adjustable cycle time between 60 and 120 seconds
- Built in countdown timer or continuous mode
- Focus Therapy • Pause Button • Compliance Meter



SC-2008-DL (E0651)

- 8 chamber compression therapy
- Unilateral or bilateral capability
- Adjustable pressure between 20 mmHg and 120 mmHg
- Adjustable cycle time between 60 and 120 seconds
- Built in countdown timer or continuous mode
- Focus Therapy • Pause Button • Compliance Meter



SC-4004-DL (E0652)

- 4 chamber compression therapy
- Unilateral or bilateral capability
- Touch Screen LCD
- Adjustable pressure between 20 mmHg & 120 mmHg in all chambers
- Adjustable cycle time between 60 and 120 seconds
- Built in countdown timer or continuous mode
- Focus Therapy • Pause Button • Compliance Meter



SC-4008-DL (E0652)

- 8 chamber compression therapy
- Unilateral or bilateral capability
- Touch Screen LCD
- Adjustable pressure between 20 mmHg & 120 mmHg in all chambers
- Adjustable cycle time between 60 and 120 seconds
- Built in countdown timer or continuous mode
- Focus Therapy • Pause Button • Compliance Meter

- 3 Year Warranty on all pumps
- Prescription required
- Contraindications apply



Bio Compression Systems
120 West Commercial Avenue
Moonachie, NJ 07074
Phone: 201-939-0716/ Fax: 201-939-4503
email: biosystems@biocompression.com
www.biocompression.com