

Bio Compression Systems' Biocomfort Garments

Made of high quality 200 Denier Nylon Oxford coated with 3 mils of Polyurethane
Available in a variety of sizes and styles including: over the shoulder,
adjustable, custom and amputee.

Upper Extremity/4-Chamber Garment Sizing Chart

Part Number	Size	Circumferences					Lengths	
		Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips
GS-3035-S	Small	23.5"	22"	19.5"	17.5"	15"	24"	26.5"
GS-3035-M	Medium	27"	26"	24"	21"	18"	28"	31"
GS-3035-L	Large	25.5"	24"	21.5"	19"	16"	31"	35.5"
GS-3035-SH/S	Small	23.5"	22"	19.5"	17.5"	15"	24"	26.5"+8"
GS-3035-SH/M	Medium	27"	26"	24"	21"	18"	28"	31"+8"
GS-3035-SH/L	Large	25.5"	24"	21.5"	19"	16"	31"	35.5"+8"

NOTE: Measurements above are in **Inches** and are the maximum allotted point on limb

Lower Extremity/4-Chamber Garment Sizing Chart

Part Number	Size		Circumferences					Lengths		
			Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease
GS-3045-H	Half	Standard	N/A		29"	27"	23.5"	22"	19"	NA
GW-3045-H		Wide			30.5"	28"	29"	23"		
GA-3045-H		Adjustable			32"	30"	23.5"	22"		
GWA-3045-H		Adj. Wide			33.5"	31"	29"	23"		
GN-3045-S	Small	Narrow	25.5"	25"	22"	20"	18"	19.5"	25.5"	27"
GS-3045-S		Standard	31"	30"	28"	24.5"	21.5"	19.5"		
GW-3045-S		Wide	35.5"	34"	30.5"	28"	25.5"	21.5"		
GA-3045-S		Adjustable	34"	33"	31"	27.5"	21.5"	19.5"		
GWA-3045-S		Adj. Wide	38.5"	37"	33.5"	31"	25.5"	21.5"		
GN-3045-M	Medium	Narrow	27.5"	25.5"	23"	20"	18"	19.5"	28"	32"
GS-3045-M		Standard	32"	30"	28"	24.5"	21.5"	19.5"		
GW-3045-M		Wide	36.5"	34"	30"	28"	25.5"	22.5"		
GXW-3045		Extra Wide	44"	41"	36"	34"	30"	22.5"	29"	32"
GA-3045-M		Adjustable	35"	34"	31"	27.5"	21.5"	19.5"	28"	32"
GWA-3045-M		Adj. Wide	39.5"	36.5"	33.5"	31"	25.5"	22.5"		
GXWA-3045-M		Adj. Extra Wide	47"	44"	39"	37"	30"	22.5"	29"	32"
GN-3045-L	Large	Narrow	28"	26"	24"	21"	18"	19.5"	32"	36"
GS-3045-L		Standard	33"	32"	30"	26"	21.5"	19.5"		
GW-3045-L		Wide	37"	36"	32"	29"	24"	20.5"		
GA-3045-L		Adjustable	37"	36"	33"	29"	21.5"	19.5"		
GWA-3045-L		Adj. Wide	40"	38"	35"	32"	24"	22.5"		

NOTE: Measurements above are in **Inches** and are the maximum allotted point on limb



Bio Compression Systems
120 West Commercial Avenue
Moonachie, NJ 07074
Phone: 800-888-0908/ Fax: 201-939-4503
email: biosystems@biocompression.com
www.biocompression.com