


PressureGuard EASY AIR®



Easy Air's patented "Air Diffusion Matrix" design maximizes removal of excess moisture.



Therapy Provided

- ✓ Microclimate Management ["Low Air Loss" (LAL)]
- ✓ Alternating Pressure
- ✓ Powered Pressure Redistribution

QUICK FACTS

Indications:

- Excessive perspiration or maceration
- Treatment of Stage 1-4 pressure injuries
- High risk patients requiring aggressive prevention
- Patients with recent flaps/grafts

Key Advantages:

- Documented superiority in moisture removal
- Won't deflate in power outage
- Includes alternating pressure mode
- PressureGuard® design: stability, safety, anti-shearing

Item #	Description
L7535-29	Easy Air System 75"L x 35"W x 7"H
* L7835-29	Easy Air System 78"L x 35"W x 7"H
L8035-29	Easy Air System 80"L x 35"W x 7"H
L8435-29	Easy Air System 84"L x 35"W x 7"H
* L8032-29	Easy Air System 80"L x 32"W x 7"H
* L8039-29	Easy Air System 80"L x 39"W x 7"H
8000	Easy Air control unit only
CLT-L7535	Easy Air Coverlet 75"L x 35"W
* CLT-L7835	Easy Air Coverlet 78"L x 35"W
CLT-L8035	Easy Air Coverlet 80"L x 35"W
CLT-L8435	Easy Air Coverlet 84"L x 35"W
* CLT-L8032	Easy Air Coverlet 80"L x 32"W
* CLT-L8039	Easy Air Coverlet 80"L x 39"W

* Allow up to 10 working days lead time for this size.

Complete System Includes Mattress, Control Unit, and Coverlet.
To order Coverlet only, add "CLT-" to part number (example: CLT-L8035).


Medicare Approved:	E0277 (Group II)
Weight limit:	500 lbs.
Warranty:	2 years, non-prorated
Mattress Weight:	22 lbs.
UL listed:	Yes

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.

Clinicians choose microclimate management to address excessive perspiration and heat that may cause maceration and lead to skin breakdown. The PressureGuard Easy Air is independently documented to outperform every leading microclimate management system (LAL) in this key performance area. The exclusive Air Diffusion Matrix™ technology and patented PressureGuard air/foam design provide microclimate management (LAL) without compromising pressure redistribution and support.

Easy Air maintains inflation in a power outage. Its stable surface won't "swallow up" the user. And its Safety Edge™ won't collapse, making transfers safer. Doesn't make the user too hot or too cold, either. Ultra-simple to set up, clean and use. Selectable comfort level. "On/Off" buttons control air flow and alternating pressure features separately.

System includes mattress, controller unit and machine washable coverlet.



Data Source:
A Proposed Method for Quantifying Low-Air-Loss Mattress Performance by Moisture Transport, Richard S. Figliola PhD, PE, Clemson University, Clemson, SC, *Ostomy/Wound Management*, 2003;49(1):32-41.

