

ORTHO8™

CRYO8 HIP CRYO HIP ROM ORTHOSIS WITH 2 GELS

PRODUCT FEATURES

- Hip and thigh Velcro Straps can be cut to size for a more customized fit
- Brace is fitted with a removable gel pack for optional cryo therapy
- Brace will control ROM (range-of-motion) at 15-degree increments
- ROM hinge and brace is also removable, allowing it to be used for support purposes

INDICATIONS

- Hip revision survey/injury to relieve pain
- Protection and support for hip instability problems
- Arthroscopic hip procedures
- Labral repair, with or without gluteous medius repair
- Post-op proximal hamstring repair
- Patients at risk to dislocate

NEW PRODUCT LAUNCH



Order P/N: 08-0001 (Left) 08-0002 (Right)
HCPC: L1686

APM006 Rev. A

www.ortho8.com

info@ortho8.com

2217 Plaza Drive
Rocklin, CA 95765

888-963-6265