



ORTHO8™

CRYO8 KNEE CRYO PNEUMATIC ROM KNEE ORTHOSIS WITH 2 GELS

PRODUCT FEATURES

- ROM hinges with malleable aluminum uprights offering full flexion abd extension control with settings at 0°, 15°, 30°, 60°, and 90°
- Anatomical knee wrap designed to combine cold therapy and pneumatic compression
- Removable pump bulb that detaches after desired compression is achieved
- Plush foam laminate with contact closure straps
- Removable gel pack with tricot lining and non-toxic gel

INDICATIONS

- ACL repair
- Post total joint surgery
- Meniscal tears
- Patella realignment
- Rehabilitation



Order P/N: 40-4001
HCPC: L1833

APM004 Rev. A

www.ortho8.com

info@ortho8.com

2217 Plaza Drive
Rocklin, CA 95765

888-963-6265