



ORTHO8™

CRYO8 SHOULDER CRYO PNEUMATIC SHOULDER ORTHOSIS WITH 2 GELS

PRODUCT FEATURES

- Anatomical shoulder wrap designed to combine cold therapy and pneumatic compression
- Removable pump bulb that detaches after desired compression is achieved
- Plush foam laminate with contact closure straps
- Removable gel pack with tricot lining and non-toxic gel
- Universal sizing fits up to a 48" chest circumference
- Additional universal extension straps available

INDICATIONS

- Post-op shoulder
- AC joint separations
- Bursitis
- Sport injuries, therapy, and athletic training
- Pain relief
- Reduce swelling and muscle spasms



Order P/N: 08-50-5001
HCPC: L3670

APM005 Rev. A

www.ortho8.com

info@ortho8.com

2217 Plaza Drive
Rocklin, CA 95765

888-963-6265