

## MedGyn Ultrasonic Pocket & Table-Top Dopplers

*The MedGyn Ultrasonic Pocket Doppler and Table-Top Doppler meet the routine examination requirements of obstetricians, midwives and pregnant women. Both dopplers are cost effective, easy to use, light weight and portable.*

### Features & Benefits

- LCD display
- Built-in audio recorder

#### Pocket Doppler

- Accurate FHR detection with clear sound
- Interchangeable 2MHz, 3MHz, 4MHz, 5MHz, 8MHz waterproof probes
- Automatic and manual FHR counting

#### Table-Top Doppler

- Convenient remote control probe
- Dual high-fidelity speakers
- Accurate FHR calculation
- Portable, lightweight design
- SD6 Model has a infrared wireless probe



Pocket Doppler



Table-Top Doppler SD5-wired



Table-Top Doppler SD6-wireless



# Specifications

	Pocket Doppler	Table-Top Doppler
Display	45mm x 25mm LCD	62mm x 44mm LCD
Heart Rate Range	50 ~ 210bpm	50 ~ 210bpm
Resolution	1bpm	1bpm
Audio Sampling Frequency	4KHz/8KHz	4KHz
Recording Length	480s/240s	240s
Audio Output Power	0.5W	2W
Size	32mm x 85mm x 138mm	265mm x 155mm x 180mm
Weight	290g (with battery)	1.7kg (with battery & probe)



Pocket Doppler with bag



Carry Bag for SD5/SD6

## Pocket Doppler

- 017102 Pocket Doppler w/ 2MHz Probe, Rechargeable
- 017127 Pocket Doppler w/ 3MHz Probe, Rechargeable
- 017101 Pocket Doppler w/ 2MHz Probe
- 017126 Pocket Doppler w/ 3MHz Probe
- 017128 2MHz probe
- 017129 3MHz probe
- 017130 4MHz probe

## Table-Top Doppler

- 017150 SD5 Table-Top Doppler w/2MHz wired probe
- 017151 SD6 Table-Top Doppler w/2MHz wireless probe
- 017155 2MHz wired probe for SD5
- 017156 2MHz wireless probe for SD6
- 017157 3MHz wired probe for SD5
- 017158 3MHz wireless probe for SD6
- 017159 Carry Bag for SD5/SD6
- 017160 Rolling Stand for SD5/SD6
- 017199 Doppler Gel 8 oz

## TRUSTED PRODUCTS FOR PREMIER WOMEN'S HEALTHCARE

### MEDGYN PRODUCTS, INC.

100 W. Industrial Rd. Addison, IL 60101 USA

**t:** +1 630.627.4105 toll-free: 800.451.9667

**f:** +1 630.627.0127

**e:** info@medgyn.com

**w:** medgyn.com

ULTRADPG-PF2

