



## GV-3010 Elite-8-Chamber Bio Vest

This 8 chambered all in one sleeve provides tissue fluid movement from the arm all the way to the thoracic duct. With a chamber located directly at the axillary area, scar tissue and axillary nodes are effectively addressed. Chest edema is not a problem with the Bio Vest. If the thoracic duct is blocked and/or not draining adequately, fluids will be redirected in a downward motion thereby creating new channels. This process, called lymphangiogenesis, is the magic performed by this sleeve and the ultimate goal for tissue fluid movement. Donning the garment is simple and the initial fitting insures a proper fit every time with no migration while in use. The Elite-8-Chamber Bio Vest provides the best possible solution for those with difficult conditions.

### Key Features:

- Axilla Compression
- Chest Compression
- Triple Action Fluid Movement:
  - Axilla
  - Thoracic Duct
  - Directional Flow
- Adjustable fit with no migration.
- Easy Donning and Doffing

### Indications for Use:

Primary/Congenital Lymphedema

- Milroys
- Praecox
- Tarda



Secondary Lymphedema:

- Lymph Node Dissection
- Lymph Node Obstructions
- Surgical Damage to the Lymph Vessels and/or nodes
- Radiation Therapy
- Trauma

Indications are not limited to the above, others may be determined by the treating physician.

Prescription Required





Supply Items  
E0657 – Vest  
E0668 – Arm Garment (S,M,L)

### Contraindications:

- Active Deep Vein Thrombosis (DVT) – Blood Clot
- Active Cancer – with the exception of palliative care
- Infection – Prior to 72 hours of antibiotic therapy
- Congestive Heart Failure (CHF)
- Any other condition wherein compression is not recommended

## Upper Extremity/8-Chamber Garment Sizing Chart

Part Number	Size	Circumferences					Lengths		
		Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips	Thorax
GV-3010-S-L	Small Left	25"	23.5"	21"	20"	19"	25.5"	29.5"	Adjustable 28" to 55"
GV-3010-S-R	Small Right	25"	23.5"	21"	20"	19"	25.5"	29.5"	Adjustable 28" to 55"
GV-3010-M-L	Medium Left	25"	23.5"	21"	20"	19"	28"	32"	Adjustable 28" to 55"
GV-3010-M-R	Medium Right	25"	23.5"	21"	20"	19"	28"	32"	Adjustable 28" to 55"
GV-3010-L-L	Large Left	25"	23.5"	21"	20"	19"	31"	35"	Adjustable 28" to 55"
GV-3010-L-R	Large Right	25"	23.5"	21"	20"	19"	31"	35"	Adjustable 28" to 55"

NOTE: Measurements above are in **Inches** and are the maximum allotted point on limb



Bio Compression Systems  
120 West Commercial Avenue  
Moonachie, NJ 07074  
Phone: 800-888-0908/ Fax: 201-939-4503  
email: biosystems@biocompression.com  
www.biocompression.com