



Best-in-Class Effective Lifting Range and Technology to Promote User Compliance

The Hoyer® HPL500 is a versatile, heavy-duty mobile patient lift that can easily handle most any patient handling task. The lift has features of a premium patient lift with a price you would expect to pay for a mid-range patient lift. Your patient will be safer and less anxious because they are at a greater distance from the lifting mechanism during transfer. This is accomplished by the HPL500's unique ELR (Effective Lifting Range) making it the preferred choice for caregivers.

EFFECTIVE LIFTING RANGE (ELR)



Patients will be safer and less anxious as the ELR keeps them at a greater distance from the lifting mechanism during critical transfers: chair-to-chair, bed-to-chair, floor-to-chair and chair-to-bath.

STANDARD POWERED BASE



Powered leg adjustment allows the caregiver to open and close the legs without bending, twisting or use of force.

OPTIONAL MANUAL BASE



The optional manual base enables caregivers to use their leg muscles and bodyweight to easily open and close the base legs.

OPTIONAL SCALE



Scale is simple to operate and has a "zeroing out" feature that allows you to weigh the patient during transfer process, without including the weight of the sling.



SAFE WORKING LOAD: 500 LB / 227 KG



HOYER SMART MONITOR

- ◊ Lift Cycle Counter
- ◊ Overload Protection
- ◊ Service Reminder
- ◊ Emergency

SPECIFICATION

Maximum Overall Length	54.4"	1382 mm
Minimum Overall Length	50.1"	1273 mm
Maximum Overall Height	83.7"	2125 mm
Minimum Overall Height	56.4"	1433 mm
Legs Open - External Width	49.5"	1257 mm
Legs Open - Internal Width	44.9"	1141 mm
Legs Closed - External Width	27.4"	696 mm
Legs Closed - Internal Width	22.4"	569 mm

WEIGHTS

Mast, Powered Base and Boom	98.9 lb	44.9 kg
Power Pack	6.1 lb	2.8 kg
Powered Base Total	105 lb	47.6 kg
Mast, Manual Base and Boom	100.8	45.7 kg
Manual Base Total	106.9 lb	48.5 kg

Specifications may vary due to manufacturing tolerances.