



CircuFlow™ 5208

Lymphedema Compression Therapy

Product Overview

The **CircuFlow™ 5208** is a lymphedema compression therapy pump used to help relieve the symptoms associated with chronic lymphatic and venous disorders in the extremities. Using intermittent gradient sequential pressure applied through easy-to-use multi-chambered garments, the **CircuFlow™ 5208** increases blood flow and provides extra cellular fluid clearances.

Features and Benefits:

- Safe:** Prevents pressure tunneling due to operator setup errors.
- Effective:** Two modes available to maximize treatment.
- Comfortable:** Able to skip chambers if necessary to ensure comfort.

Device Specifications:

Weight:	12 lbs.
Inflation Cycle:	Adjustable 2.5 - 30 seconds per chamber
Deflation Cycle:	12 seconds
Dimensions:	11.5" (W) x 10" (D) x 4" (H)
Cycle Time:	Adjustable 20 - 240 seconds
End Therapy Vacuum Time:	60 seconds
Therapy Time:	Adjustable 15 - 180 minutes
Electrical:	120VAC, 60Hz, 100VA MAX
Fuse Rated:	250VAC, 1.0AMP, SLO-BLO
Operation Mode:	Continuous Operation with Intermittent Loading



Item Code: P-5208
LB35.0002 Rev E 20160115