


FUSION^{2K}[™]

Featuring SelectAir[®] Therapy

Wicking Action: Air escapes through tiny holes on the top of the mattress and circulates under a vapor-permeable cover to wick moisture away from patient.



iSense Processing[™]

Advanced microprocessors continually monitor patient's immersion.

FUSION^{2K} Features

'True' Low Air Loss Technology

Allows air to escape from the top of the mattress and directly interact with moisture at the patient/bed boundary.

Pulsation (alternating pressure)

Provides two cycle settings that continually readjust the system to ensure patient comfort.

Firm Perimeter Support

For ease of transfer and safety.

User-Friendly Control Unit

Rapid inflation, easy set-up and a motor so quiet, the patient will hardly know it's there!

Zoned, "true" low air loss used in acute care hospitals nationwide with firm perimeter support for ease of transfer and safety

Weight Limit: 300 lbs. (136 kg)

Mattress Size: 80" (or 84") x 36" x 6"

Cover: Water resistant, vapor-permeable, flame resistant, biocompatible and quilted

Warranty: One Year Limited Warranty on control unit and soft goods

Medicare Code: E0277 Group 2 Support Surface

LIT-0007 | R100121

Manufactured in the USA by
MOXI[™]
ENTERPRISES
Innovative health solutions

1939 S. Vandeventer Ave.
St. Louis, MO 63110
T 1-800-942-7678
E info@moxiusa.com
W moxiusa.com