



3 Wheel Scooter

MOB1025

We offer you a mobility 3 wheels scooter that is easy to transport and easy to maneuver, without sacrificing quality. Whether you or your loved ones need help getting around the house and shopping or want the freedom to travel. This 3-wheel mobility scooter can easily maneuver through hallways and doors and is compact enough to be used while traveling.



Specifications

- What's Included
 - Vive 3-Wheeled Mobility Scooter
 - Off-board battery charger
- Colors:
 - Blue
 - Red
 - Silver
- Sizes:
 - Front Tires: 7.5"
 - Rear Tires: 7.5"
 - Overall Length: 37"
 - Overall Width of Base: 19.67"
 - Seat-to-Ground Height Range: 18.5" to 21.5"
 - Seat Dimensions: 17" x 14.5"
 - Seat-to-Deck Height Range: 16.25" - 17.25"

Distributed by



8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com



Comfortable Seating

Swivel seat and contoured backrest are height adjustable and covered with a supple black vinyl for greater comfort and easy cleaning



More Standard Features

Each scooter includes a sturdy front tiller basket, a high-intensity LED headlight and horn for added safety



Variable Speed Scooter

An intuitive control panel with the speed adjustment knob provides greater safety and control when navigating indoor and outdoor surfaces



Stable Maneuverability

Easily navigate over smooth surfaces and up and down inclines of nine degrees or less

Distributed by



8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

☎ (239) 220-5772 🌐 vivehealth.com ✉ service@vivehealth.com
📍 8955 Fontana Del Sol Way, Naples, FL 34109, USA

SR#: 2021-01-R1