

4 Wheel Mobility Scooter

MOB1027

Move Smooth Surfaces Safely - Easy to handle on smooth surfaces, indoors and out. You can move around comfortably with the Vive four-wheel mobility scooter. Designed to safely support up to 265 pounds, the mobility scooter features 8-inch tires, a top speed of 3.7 mph, and a driving distance of up to 12.4 miles per charge.



Specifications

- What's Included
 - Vive 4-Wheeled Mobility Scooter
 - Off-board battery charger
- Colors
 - Blue
 - Silver
 - Red
- Sizes:
 - Overall Length: 40"
 - Overall Width of Base: 19.68"
 - Seat-to-Ground Height Range: 18.5"-21.5"
 - Seat Dimensions: 17" x 14.5"



Variable Speed Scooter

Select the maximum speed desired with the speed adjustment dial, providing greater safety and control



More Standard Features

Each scooter includes a removable front tiller basket, a high-intensity LED headlight and horn for added safety



Comfort Swivel Seating

Easy to board and exit, the swivel seat & contoured backrest are height adjustable & covered with a supple black vinyl for comfort & easy cleaning



Easy Maneuverability

Resting on a stable base of four, flat-free, tires to easily navigate over smooth surfaces and up and down inclines of six degrees or less