



vive health innovation

vive® Compact Power Wheelchair

MOB1029S

Providing comfortable, yet stylish transportation, the Vive Compact Power Wheelchair allows you to easily navigate on surfaces safely. The durable aluminium frame is compact when folded to easily fit in standard trunk sizes. Additionally, the chair is ideal for airline and cruise ship travel.

Safely supporting up to 220 pounds, the electric wheelchair includes powerful batteries that handle mild inclines and flat surfaces up to 12 miles per charge. With a maximum speed of 4 mph, the wheelchair is easy to operate with a 360-degree joystick controller that can be positioned on the left or right height-adjustable armrest.

Specifications

- What's Included
 - 1x Compact Power Wheelchair
 - 2x Batteries
 - Electrical components
- Colors: Silver/black
- Weight Limit: 220 lbs
- Product Weight (with battery): 95 lbs
- Overall Size: 39.5" x 26" x 38"
- Top Speed: 4 mph
- Max Drive Range Per Charge: 12 miles
- Latex Free

Materials:
Led acid battery, aluminium frame



Distributed by



8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com



Easy to Maneuver

Features a repositionable, 360-degree joystick for easy maneuverability.

Comfort Seating

Secured with an onboard lap belt, the padded seat is comfortably wide and covered with a soft, breathable material.

Compact & Foldable

Remove the front leg rests and fold for easy, compact storage and travel.

Powerful Dual Motors

Smoothly navigate both flat surfaces and mild inclines.



Adjustable Foot Rests

Rotating, foldable foot rests make it easy to get in and out of the wheelchair.

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

☎ (239) 220-5772



vivehealth.com



service@vivehealth.com



8955 Fontana Del Sol Way, Naples, FL 34109, USA

SR20-5R3