

## MedGyn AL-102 & AL-401 Colposcopes

*MedGyn's AL-102 (1-step with 2 interchangeable eyepieces) & AL-401 (3-step magnification) colposcopes are easy to operate, with excellent optics and the preferred option for smaller practices looking to benefit from lost cost options. MedGyn's basic stereoscopic colposcopes have become the industry standard by combining high-quality lens systems with easy maneuverability.*

### Features & Benefits

- Easy-to-use, stable, light-weight tilt stand
- Electronically-controlled illumination
- Green filter
- Smooth action controls
- Interchangeable eyepieces

#### AL-102 Colposcope

- Interchangeable eyepieces at 10x, 13.5x, 18x magnification
- 3-position illumination settings
- Complete video and camera system adaptability

#### AL-401 Colposcope

- 3-step magnification: 5x, 10x, 20x
- Precision computer-designed 10x optics
- 300mm focal length



MedGyn AL-102



MedGyn AL-401



## Simplicity

	AL-102	AL-401
<b>Eyepiece</b>	10x, 13.5x, 18x	10x
<b>Magnification</b>	8x, 10.8x, 14.4x	5x, 10x, 20x
<b>Working Distance</b>	275mm	300mm
<b>Field of View</b>	10x = 20mm 13.5x = 15mm 18x = 12.5mm	20mm
<b>Weight</b>	28lbs, 11kg	18lbs, 8kg

## Illumination

	AL-102	AL-401
<b>Lamp</b>	6v, 2A, 12w	12V, 20 LED Light Source
<b>Main Power</b>	110V or 220V	110V/220V e
<b>Illumination Field</b>	25mm	50mm
<b>Filter</b>	Green	Green

### Colposcope

010200 MedGyn AL-102 Colposcope

010198 MedGyn AL-401 Colposcope

### Accessories & Parts

010204 AL-102 10x eyepiece (set of 2)

010205 AL-102 18x eyepiece (set of 2)

010206 AL-102 13.5x eyepiece (set of 2)

010001AL-102 Video System with video pickup and color monitor

010209 AL-102 Replacement Bulb

010207 AL-102 Replacement Fuse

## TRUSTED PRODUCTS FOR PREMIER WOMEN'S HEALTHCARE

### MEDGYN PRODUCTS, INC.

100 W. Industrial Rd. Addison, IL 60101 USA

**t:** +1 630.627.4105 toll-free: 800.451.9667

**f:** +1 630.627.0127

**e:** info@medgyn.com

**w:** medgyn.com

AL114COP-PF2

