



# MedGyn

Innovative Products for  
Ob/Gyn Professionals

## MedGyn AL-106HD™ AL-106HD-SA Digital Video Colposcopes & Image and Data Management Software

### Simplicity

*Different from traditional optical colposcopes, the AL-106HD™ and AL-106HD-SA are digital colposcopes with an LED light source to ensure ultra-bright visualization and a long life cycle. Their unique optics provide multiple filters to ensure high definition and quality images and improved visualization of blood vessels.*



AL-106HD



AL-106HD Rolling Base



AL -106HD Swing Arm



## Features & Benefits

	AL-106HD	AL-106HDSA
Pixels	1.43M	1.43M
Magnification	40x	40x
Connection	HDMI	HDMI
Light Source	White / Blue LED	White / Blue LED
Horizontal Resolution	1080	1080
Camera	CMOS	CMOS

## Technical Specifications

	AL-106HD	AL-106HD-SA
Filters	3x Green	3x Green
Working Distance 200mm - 300mm	√	√
Depth of View 5mm - 120mm	√	√
Video Capability	√	√
Field of View 60mm - 6mm	√	√
Auto/Manual Focus Function	√	√
Image Freeze Function	√	√
Main Power 110/240VC 50/60Hz	√	√
Battery Backup Option		
Wireless Signal Transmission		

## Optional Features

	AL-106HD	AL-106HDSA
Swing Arm		√
Image & Data Management Software	√	√
Desktop/ Laptop Computer System	√	√
Flat-Screen Color Monitor	HDMI Port	HDMI Port






# MedGyn Image & Data Management Software with Forensic Capability

## Simplicity

*Designed for use with MedGyn AL-106HD™ and AI-106HD™-SA Digital Video Colposcopes, the Image & Data Management Software allows for the capture and management of colposcopy images into a self-contained patient record on an easy, user-friendly and secure system.*

## Features & Benefits

- Complete system security, internal password protection
- Easily customizable to specific clinical requirements and preferences
- Includes a large offering of standard examination and diagnosis templates
- Full video capture capability
- Statistical analysis
- Image measurement & annotation
- System automatically creates variety of report options
- Software includes: capture card, security key and foot switch
- Compatible with Windows 7, Windows 8 & Windows 10



## AL-106HD Colposcope

010466SA-HD AL-106 High Definition Digital Video Colposcope with Swing Arm

010466-HD AL-106 High Definition Digital Video Colposcope

## Image and Data Management Software

010475 Image and Data Management Software (with Desktop/ Laptop and Monitor)

010476 Image and Data Management Software installed on Desktop/ Laptop (without Monitor)

010477 Image and Data Management Software

## TRUSTED PRODUCTS FOR PREMIER WOMEN'S HEALTHCARE

### MEDGYN PRODUCTS, INC.

100 W. Industrial Rd. Addison, IL 60101 USA

**t:** +1 630.627.4105 toll-free: 800.451.9667

**f:** +1 630.627.0127

**e:** info@medgyn.com

**w:** medgyn.com

AL106HDS-PF2

