

MedGyn AL-106™ and AL-106-SA™ Digital Video Colposcopes & Image and Data Management Software

Simplicity

Different from traditional optical colposcopes, the AL-106™ and AL-106-SA are digital colposcopes with an LED light source to ensure ultra-bright visualization and a long life cycle. Their unique optics provide multiple filters to ensure high quality images and improved visualization of blood vessels.



AL-106



AL-106 Swing Arm



Features & Benefits

	AL-106	AL-106-SA
Pixels	440K	440K
Magnification	40x	40x
Connection	AV	AV
Light Source	White LED	White LED
Horizontal Resolution	550	550
Camera	CCD	CCD

Technical Specifications

	AL-106	AL-106-SA
Filters	3x Green	3x Green
Working Distance 200mm - 300mm	√	√
Depth of View 5mm - 120mm	√	√
Video Capability	√	√
Field of View 60mm - 6mm	√	√
Auto/Manual Focus Function	√	√
Image Freeze Function	√	√
Main Power 110/240VC 50/60HZ	√	√
Battery Backup Option		
Wireless Signal Transmission		

Optional Features

	AL-106	AL-106-SA
Swing Arm		√
Image & Data Management Software		
Desktop/ Laptop Computer System	√	√
Flat-Screen Color Monitor	LCD Video	LCD Video




MedGyn Image & Data Management Software with Forensic Capability

Simplicity

Designed for use with MedGyn AL-106™ and AL-106-SA Digital Video Colposcopes, the Image & Data Management Software allows for the capture and management of colposcopy images into a self-contained patient record on an easy, user-friendly and secure system.

Features & Benefits

- Complete system security, internal password protection
- Easily customizable to specific clinical requirements and preferences
- Includes a large offering of standard examination and diagnosis templates
- Full video capture capability
- Statistical analysis
- Image measurement & annotation
- System automatically creates variety of report options
- Software includes: capture card, security key and foot switch
- Compatible with Windows 7, Windows 8 & Windows 10



AL-106 Colposcope

- 010466 AL-106 Digital Video Colposcope
- 010466SA AL-106 Digital Video Colposcope with Swing Arm

Image and Data Management Software

- 010475 Image and Data Management Software (with Desktop/ Laptop and Monitor)
- 010476 Image and Data Management Software installed on Desktop/ Laptop (without Monitor)
- 010477 Image and Data Management Software

TRUSTED PRODUCTS FOR PREMIER WOMEN'S HEALTHCARE

MEDGYN PRODUCTS, INC.

100 W. Industrial Rd. Addison, IL 60101 USA

t: +1 630.627.4105 toll-free: 800.451.9667

f: +1 630.627.0127

e: info@medgyn.com

w: medgyn.com

AL106SAM-PF2

