

# LUMEX®

A Graham-Field® Family Brand

# Three Position Recliners

- The Standard of the Industry and the first recliner specifically designed to improve the comfort and mobility of residents in extended care facilities. The Lumex Three Position Recliner easily adjusts into any one of three positions: sitting, TV or full recline. Recliners offer numerous therapeutic benefits including position changes that can aid in circulation and ease discomfort from respiratory, heart and other ailments.
- A position lock allows only the caregiver to secure the chair in any one of the three easily attained positions. Recliner encourages natural relaxation and improved blood circulation.
- Maximum Weight Capacity: 250 lb EVENLY DISTRIBUTED

### Warranty

Welded Frame: Lifetime Warranty (Defined as Ten (10) Years)

Mechanical Components: Three (3) Years

Upholstered Components: One (1) Year


- The molded activity table can be fixed in position or left unsecured for easy exit from the recliner.
- Ergonomically designed seat and back with varied density foam construction provides correct support. Articulating headrest quickly adjusts to maximize individual comfort requirements.
- For safer, easier transport, the Lumex Three Position Recliner features two, four inch swivel casters with wheel locks and five inch fixed front wheels.
- Durable welded steel epoxy coated frame is easy to clean and maintain.

Casters: One (1) Year


\*Subject to Warranty Statement Limitations



Available Colors*			UOM
574G401	FR574G401	Vintage Gold	1 ea
574G409	FR574G409	Warm Taupe	1 ea
574G427	FR574G427	Blue Ridge	1 ea
574G432	FR574G432	Imperial Blue	1 ea
574G454	FR574G454	Royal Blue	1 ea
574G851	FR574G851	Doe Skin	1 ea
574G857	FR574G857	Jade	1 ea
574G863	FR574G863	Rosewood	1 ea



\* FR models are Fire Rated Colors (Meets California Technical Bulletin 133 Flammability Standards)



Locked

Unlocked



www.grahamfield.com

© 2018 GF Health Products, Inc.  
Atlanta, GA All rights reserved.  
GF0900049RevB18