


# Protekt™ 200

Pressure Redistribution Foam Mattress



 **Proactive**  
medical products

Medicare  
Code  
E0184

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

[www.proactivemedical.com](http://www.proactivemedical.com) • [sales@proactivemedical.com](mailto:sales@proactivemedical.com)

The Protekt™ 200 is an economical multi-tiered pressure redistribution mattress. This mattress is ideal for moderate to high risk patients.


## FEATURES AND BENEFITS:

- Multi-tiered pressure redistribution mattress with top layer of die-cut high density foam in the head and torso section provides superior pressure redistribution.
- Bottom layer of high density foam provides a firm support base.
- Heel section is multi-tiered, the top layer consist of individually responsive cells for superior pressure redistribution.
- PTSS Rail System provides additional support, allowing easy patient transfer and edge of bed sitting.
- Top cover is a 4-way stretch nylon which reduces perspiration and provides a low shear friction surface that is anti-microbial and water resistant.
- Designated head and foot sections.
- Anti-Contamination Flap covers the zipper reducing the possibility of fluids contaminating the foam core.
- Fire sleeve to meet Federal Fire Code 16 CFR part 1633.
- Weight Capacity: 300 Lbs.
- Warranty: 5 year, Non Prorated Warranty.

81021	Protekt 200	Pressure Relieving Foam Mattress	36"x76"x6"
81022	Protekt 200	Pressure Relieving Foam Mattress	36"x80"x6"
81023	Protekt 200	Pressure Relieving Foam Mattress	36"x84"x6"
81024	Protekt 200 Plus	Pressure Relieving Foam Mattress w/3" Raised Rails	36"x76"x9"
81025	Protekt 200 Plus	Pressure Relieving Foam Mattress w/3" Raised Rails	36"x80"x9"
81026	Protekt 200 Plus	Pressure Relieving Foam Mattress w/3" Raised Rails	36"x84"x9"



**Optional raised side rails**



270 Washington Street Mount Vernon, NY 10553

**Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005**

[www.proactivemedical.com](http://www.proactivemedical.com) • [sales@proactivemedical.com](mailto:sales@proactivemedical.com)