

Short-Wave™

Wheelchair Seat and Back Cushion

- Designed for standard wheelchairs
- Available in 16", 18" or 20 widths
- Readily adaptable to any seating application



Reduces Shear

Unique surface design minimizes the effects of shearing.

Assists Positioning

Natural lumbar and thigh supports help to position users comfortably.

User-Friendly

One-piece construction, lightweight and durable. Wipe-clean, non-skid cover attaches easily to the chair.



Medicare approved: Bills under **BOTH E2601 (Seat) and E2611 (Back)** codes concurrently.

SPANAmerica
Innovative Solutions.

Short-Wave™

Simple Secure Attachment

Two elastic hook and loop straps located at the back and seat help anchor the cushion to the chair.



Fluid-repellant Cover

Removable zippered cover features non-slip top fabric and durable non-skid bottom. Wipes clean with neutral suds or approved cleaning agent.

Proprietary Wave Profile



Unique segmented fins are specifically engineered to address the effects of shearing that are typical of institutional seating.

Patented design has two distinct zoned orientations to reduce shear in multiple directions for specific seating areas. (Use of cushion also helps to reduce pressure.)

Item	Description (cushion width)	Qty
SWAVE16-04	16" Short-Wave Cushion w/Cover	4/cs
SWAVE16-01	16" Short-Wave Cushion w/Cover	1/cs
SWAVE-04	18" Short-Wave Cushion w/Cover	4/cs
SWAVE-01	18" Short-Wave Cushion w/Cover	1/cs
SWAVE20-04	20" Short-Wave Cushion w/Cover	4/cs
SWAVE20-01	20" Short-Wave Cushion w/Cover	1/cs

Warranty: 18 Months

Weight Limit: 300 Pounds

Designed to fit 16" or 18" deep seat

SPANAmerica
Innovative Solutions.