


Prevent Wandering

Anti-Wandering Door Systems

New & Affordable!

For Homes & Facilities of any size!



Affordable – No Wiring – No Electricians

MFR#	PRODUCT DESCRIPTION
TL-4005SYS	Single Door Bar Monitor System (1 year warranty)
TL-4012	Resident Wristband Transmitter (6 month warranty)
TL-4014W	Wall Mount Reset Button (1 year warranty)
AC-04	AC Power Adapter


Optional Additional Components:



TL-2025
Magnetic
Door Lock



TL-4015 Central
Monitor Displays
Which Resident/Device
Triggered Alert



All-in-One

This "All in One" Anti-Wandering Door Bar System functions with just a door alert and a resident wristband. Simply mount the door bar by any exit or doorway in your care facility, plug it in, and place wristband on resident. When a resident wearing a wristband gets within 8 to 14 feet of the door bar or attempts to go through the doorway, the door bar will alert audibly and visually.

Simple Installation

1. Mount Door Bar
2. Activate Wristband Transmitter
3. Test

Call Us Today! (800) 650-3637 | www.smartcaregiver.com

