

Wireless Central Monitoring

Build your own perfect system!

Starter Kit: 433-PRB-SYS

- 433-CMU-40
- 433-NC (3x)
- 433-PRB
- AC-Adaptor



433-PRB
Caregiver Pager
w/ Reset Button

433-NC (3)
Nurse Call Buttons


Add-Ons:



Option 1



Option 2



Option 3

Option 1

Chair Pad

GCT-WI 10"x 15" Weight Sensing Chair Pad

Option 2

Bed Pad

GBT-RI 10"x 30" Weight Sensing Bed Pad

Option 3

Floor Mat

FMT-07C 24"x 48" Weight Sensing Floor Mat

Central Monitor Key Features:

- Alarms audibly & visually at central location and portable caregiver pager(s).
- Compatible with all 433 Wireless Components.
- 40 & 60 Channel Capacity Models available.
- Range up to 300 ft.
- Systems may be pre-programmed prior to delivery or easily re-programmed using the provided instructions.
- Batteries not included (4-C)
- AC adaptor included (AC-04)



Call Us Today! (800) 650-3637 | www.smartcaregiver.com